### Safe Personal Effective



# Use of custom silicone putty for nasal tube stabilisation in Orthognathic surgery.

Mr G.Chiu (Consultant OMFS), E.Morley (DCT2 OMFS) Department of Oral & Maxillofacial Surgery

East Lancashire NHS Trust







Fig. 1. Anchoring of the tube on the forehead with the silicone mould and tape after nasotracheal intubation. The putty used for this impression was Express<sup>TM</sup> VPS impression material (3M ESPE, MN, USA), and the nasotracheal tube was a Parker Flex-Tip<sup>TM</sup> cuffed nasotracheal tube (Parker Medical, CO, USA).

# Background

- Several factors contribute to the successful outcomes of orthognathic surgery, including the undistorted position of the nose during surgery which is the main focus in this presentation.
- Nasal intubation is predominantly used in airway management in orthognathic surgery.
- Where surgical manipulation of the mid/upper face is required, nasal intubation may interfere with accuracy of surgical treatment as it can uncontrollably distort the nose.
- Submento-tracheal intubation may therefore be a preferred method by some surgeons, however an unnecessary facial scar results .
- A method of Anchoring of a nasotracheal tube using flexible silicone putty was described in the BJOMS in 2016 as a method to reduce intraoperative shifting of the nasotracheal tube resulting in necrosis and ulceration of the nasal alar skin.
- For orthognathic surgery in RBH we have adopted this technique in attempt to reduce uncontrolled alterations of nasal proportions.

## Aims

- To assess an alternative method involving use of silicone dental impression putty for tube stabilisation
- To assess if there is a reduction in uncontrolled alterations of naso-labial angle and inter-alar width, compared to standard nasal intubation.







### Method

- Comparisons of pre-operative and post-operative inter-alar width and naso-labial angle were made between 2 groups of patients.
- Both groups underwent nasal intubation.
- One group of patients had securing of the nasal tube to the forehead with conventional methods of swabs & tape.
- In the second group, the custom silicone-putty device was used to secure & stabilise the nasal tube.
- Measurements of naso-labial angles and inter-alar widths were taken pre-intubation and again post-intubation to determine a difference.
- The differences were then compared between the two groups to assess potential benefits of the additional method of stabilisation.



### Mann Whitney U Test

- The Mann Whitney U Test was used to assess for significant difference in the data collected.
- Null hypothesis established.
- Significance value was set at 0.05 (5%)
- USTAT value calculated for both data sets
- Lowest USTAT value of the 2 data sets taken
- Critical value determined at 0.05 significance '41' (UCRIT – using Mann Whitney U chart )
- Significant difference →demonstrated when U<sub>STAT</sub> </= U<sub>CRIT</sub> (41)
- Null hypothesis rejected if significant difference is shown.

#### CRITICAL VALUE CHART

(SIGNIFICANCE VALUE AT 5%)

n1\ <sup>n2</sup>	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
2							0	0	0	0	1	1	1	1	1	2	2	2	2
3				0	1	1	2	2	3	3	4	4	5	5	6	6	7	7	8
4			0	1	2	3	4	4	5	6	7	8	9	10	11	11	12	13	14
5		0	1	2	3	5	6	7	8	9	11	12	13	14	15	17	18	19	20
6		1	2	3	5	6	7	10	11	13	14	16	17	19	21	22	24	25	27
7		1	3	5	6	8	10	12	14	16	18	20	22	24	26	28	30	32	34
8	0	2	4	6	7	10	13	15	17	19	22	24	26	29	31	34	36	38	41
9	0	2	4	7	10	12	15	17	20	23	26	28	31	34	37	39	42	45	48
10	0	3	5	8	11	14	17	20	23	26	29	33	36	39	42	45	48	52	55
11	0	3	6	9	13	16	19	23	26	30	33	37	40	44	47	51	55	58	62
12	1	4	7	11	14	18	22	26	29	33	37	41	45	49	53	57	61	65	69
13	1	4	8	12	16	20	24	28	33	37	41	45	50	54	59	63	67	72	76
14	1	5	9	13	17	22	26	31	36	40	45	50	55	59	64	67	74	78	83
15	1	5	10	14	19	24	29	34	39	44	49	54	59	64	70	75	80	85	90
16	1	6	11	15	21	26	31	37	42	47	53	59	64	70	75	81	86	92	98
17	2	6	11	17	22	28	34	39	45	51	57	63	67	75	81	87	93	99	105
18	2	7	12	18	24	30	36	42	48	55	61	67	74	80	86	93	99	106	112
19	2	7	13	19	25	32	38	45	52	58	65	72	78	85	92	99	106	113	119
20	2	8	13	20	27	34	41	48	55	62	69	76	83	90	98	105	112	119	127

## Inter-alar Width Data

- Average difference IAW with no putty  $\rightarrow$  2.38mm
- Average difference IAW with putty  $\rightarrow$  1.35mm
- Use of putty reduces alterations of IAW's by average of  $\sim 1.03$ mm
- Mann Whitney U Test  $\rightarrow$ 
  - USTAT 42.5 > UCRIT 41
  - Significant difference <u>NOT</u> demonstrated for IAW

<u>Putty</u> Inter-alar Width						
Pre-Intubation (mm)	Post-Intubation (mm)	Difference (mm)				
29	30	1				
32	33	1				
31	34	3				
35	35	0				
30	31	1				
30	34	4				
33	33	0				
35	35	0				
33	33	0				
34	34	0				
33	35	2				
36	35	2				
32.5	36	3.5				

<u>No Putty</u> Inter-alar Width							
Pre-Intubation (mm)	Post-Intubation (mm)	Difference (mm)					
36	38	2					
36.5	37	0.5					
30	33	3					
33	36	3					
35	38	3					
32	34	2					
35	40	5					
34	38	4					
34	34	0					
28	34	6					
36	38	2					
42	40	2					

# Naso-labial Angle Data

- Average difference NLA with no putty  $\rightarrow$  19.17 degrees
- Average difference NLA with putty  $\rightarrow 0.85$  degrees
- Use of putty reduces alterations of NLA's by average of  $\sim 18.32$  degrees
- Mann Whitney U Test  $\rightarrow$ 
  - USTAT 0 < UCRIT 41
  - Significant difference <u>IS</u> demonstrated for NLA

<u>Putty</u> Naso-labial angle							
Pre-Intubation (degrees)	Post-Intubation (degrees)	Difference (degrees)					
99	95	4					
96	96	0					
95	96	1					
95	95	0					
88	95	7					
115	120	5					
89	87	2					
90	89	1					
115	115	0					
98	98	0					
88	88	0					
109	108	1					
94	100	6					

<u>No Putty</u> Naso-labial angle							
Pre-Intubation (degrees)	Post-Intubation (degrees)	Difference (degrees)					
99	120	21					
121	138	17					
112	125	13					
115	124	9					
115	138	23					
111	145	34					
100	135	35					
110	120	10					
101	115	14					
99	125	26					
98	115	17					
112	123	11					

# Summary



- Reduced alterations in both IAW and NLA measurement was shown on average with use of putty.
- Significant difference was <u>NOT</u> demonstrated in IAW measurements with use of putty stabilisation.
- Significant difference <u>WAS</u> demonstrated in NLA measurements with use of putty stabilisation.
- Silicone dental impression putty is relatively cheap to use.
- Can be effectively used to prevent nasal septa distortion during orthognathic surgery.
- Aids in the maintenance of facial symmetry and bony contours, avoiding the need to use submental intubation.

# References

- 1) Ogawa T, et al., Anchoring of a nasotracheal tube using flexible silicone putty. Letters to the Editor / British Journal of Oral and Maxillofacial Surgery 54 (2016) 475–480
- Beck J, Johnston K. Anaesthesia for cosmetic and functional maxillofacial surgery. Continuing Education in Anaesthesia Critical Care & Pain, Volume 14, Issue 1, February 2014, Pages 38 -42, https://doi.org/10.1093/bjaceaccp/mkt027
- 3) Mitchell V. Anaesthesia for orthognathic surgery | Oxford Textbook of Anaesthesia for Oral and Maxillofacial Surgery | Oxford Academic (oup.com), pages 240-252.June 2010https://doi.org/10.1093/med/9780199564217.003.0018
- 4) Malhotra N, Retromolar intubation a simple alternative to submental intubation - Malhotra – May 2006 - Anaesthesia - Wiley Online Library 61(5) pages 515-516
- 5) Altman JI, Oeltjen JC. Nasal deformities associated with orthognathic surgery: analysis, prevention, and correction. J Craniofac Surg. 2007 Jul;18(4):734-9. doi: 10.1097/SCS.0b013e3180684328. PMID: 17667658.
- 6) O'Loughlin E, How To... Perform the Mann-Whitney U Test (By Hand) -YouTube